

Pure LandingPad™ Endoscope Protective Surface Liner


Place on stainless steel surfaces where endoscopes are staged for pre-cleaning

Made from ¾ inch thick soft foam blend that resists hospital-grade chemicals and protects expensive devices and parts

Won't absorb liquids

Perforations facilitate the removal of liquids and debris from work surfaces

18X24 size




Pure LandingPad™ Endoscope Protective Sink Liner

Place at bottom of sink to protect delicate instruments

Non-toxic, eco-friendly, anti-static material

Medical grade silicone

20x20 size



PROTECTING
Precious Instruments

pure
PROCESSING

TO SEE OUR VIDEO DEMONSTRATION GO TO
www.pure-processing.com/video/

For more information on the Pure Station™
Surgical Devices Pre-cleaning System, call: **877.718.6868**