

CLEAR · SIMPLE · CONSISTENT

PureChannel™ Flush Instrument Channel Cleaner and Foam

Every aspect of this concentrated formulation is designed to help make pre-cleaning expensive medical devices easy, effective, consistent and economical. The triple-enzyme formula is specifically designed to help remove the typical organic soils that stick to hard-to-clean internal channels and external surfaces after patient use, whether it's used as a diluted liquid cleaner or a concentrated foam application.



- Double-duty formula (liquid and foam) eliminates the need for two different products
- Single dilution rate of 1/4 oz. per gallon helps eliminate confusion, inconsistency and waste
- Concentrated foam* breaks down soils, maintains moisture and facilitates cleaning
- Regular use can help reduce the need for prolonged and intensive manual cleaning
- Thoroughly tested and determined to be equally effective in a wide range of water conditions
- Available in two ergonomic sizes: a one-gallon jug with handle, and a lightweight 32-ounce bottle with a one-hand neck design for controlled dilution and easy spraying

**Optional spray nozzle available separately.*

pure
PROCESSING

For more information on FlexiPump,™ PureStation,™ and PureSteel™ ergonomic cleaning systems and products, please call: 877.718.6868 or, visit: www.pure-processing.com