

# PureStation™ Dedicated Eye Sink

## Cleaning Area & Cleaning Equipment All In One

The PureStation eye sink serves as a designated cleaning area for channeled eye surgery devices. The pump can flush up to three devices at once while pushing up to 1000cc every minute. A chemical mixing well and thermometer strip ensure that enzymatic detergent is diluted and maintained at the proper temperature and amount. The dedicated eye sink is permanently labeled with specific instructions for eye instruments, and can be used on a countertop, in a deep sink or on a mobile cart as an additional sink. Solution and rinse containers allow for the proper amount of fluids for each cleaning, and an accessory container keeps cleaning brushes and other accessories organized and secure. In addition, the sink is lightweight and easy to carry, with side handles that also serve as overflow prevention outlets.



ANSI/AAMI ST 79: 2010/Annex N states that a designated cleaning area and equipment dedicated to the cleaning of intraocular surgical instruments should be identified. Intraocular surgical instruments should be processed separately from general surgical instruments and equipment, to reduce the potential for cross contamination by material or residue from general surgical instruments.