

PureStation™ Tiered Carts

For device soaking, pre-cleaning, transport, and vertical storage

These multi-purpose tools can serve as transport vehicles, to collect soiled minimally invasive surgical, robotic, and endoscopic devices from the point of use. Once they reach a reprocessing area, they can be transformed in multiple ways, into:

- A staging area for pre-soaking two, three or four levels of soiled devices at a time
- A two-step reprocessing station, soaking below, then moving up to pre-cleaning and flushing in the top-tier PureStation™ Sink Insert with an optional FlexiPump™ Independent Flushing System
- A dedicated staging and reprocessing station for specific devices
- A safe, covered temporary storage location when space is tight or reprocessing volume is at peak levels.

Each sink insert can be fitted with a FlexiPump Independent Flushing system, for more automated flushing capacity

Inserts can also be removed and placed in a sink for ergonomic reprocessing

Optional stainless steel top adds counter space to a department if needed

Each sink insert has built-in overflow handles, temperature monitors for solutions, volume level indicators, and drain catches to retain small parts

Tiered carts and PureStation Sink Inserts can be cleaned and disinfected in cart washers

Each cart comes with the appropriate number of sink inserts and secure tight-fitting biohazard lids

Comes in two, three or four tier configurations

