

**APPLICATION**

PureSteel™ Task Lights help staff visually inspect the external surfaces of instruments for damage, and irregularities, while also providing targeted illumination. Includes multiple installation options.

**DESCRIPTION**

PureSteel™ Task Lights provide staff with greater illumination in reprocessing areas such as decontamination and prep and pack.

Features an on/off switch, and dimmable lighting. Both styles of light feature a 27" gooseneck arm for precise, controlled magnification. Water resistant to work in moist environments. Choose from a LED light source, or the halogen style with a 2x magnifier.



PureSteel Halogen Magnifier Task Light



PureSteel LED Task Light



PureSteel Desktop Halogen Magnifier Task Light

**The Selections Checked Below Apply to This Equipment**

**PURESTEEL™ TASK LIGHTS**

- PureSteel LED Task Light
- PureSteel Halogen Magnifier Task Light
- PureSteel Desktop Halogen Magnifier Task Light
- Mounting Bracket

Facility

Contact Name

Representative Name

Contact E-mail

Representative E-mail