

PureSteel™ Mobile Soaking Station

The PureSteel™ Mobile Soaking Station creates a dedicated space for soaking instruments during peak reprocessing hours. The Mobile Soaking Station easily accommodates robotic instruments to horizontally soak, and is built on casters for transport and mobile soaking needs.

The PureSteel™ Mobile Soaking Station Organizational Pegboard adds greater organization and installation options for the Mobile Soaking Station.

- Economical solution for 30-minute robotic soaking per IFUs
- Create extra sink spaces for reprocessing, or establish dedicated areas for specific instruments
- Create dedicated spaces for soaking during peak reprocessing time
- Easily accommodates robotic instruments
- Large enough to accommodate orthopedic sets
- Small footprint conserves space
- Casters allow for mobility



pure
PROCESSING

For more information about PureClear™ Inspection
Scopes and the Visualization Station:

VISIT www.pure-processing.com OR CALL 877.718.6868