

Mobile Sink



pure
PROCESSING

877.718.6868

pure-processing.com

October 14th, 2015

Mobile Sink



OSHA Central Supply, “use carts with large, low rolling, low resistance wheels, that can easily roll over mixed flooring as well as gaps between elevators and hallways.”

Our mobile cart creates another work station and is available with or without FlexiPump™ Independent Flushing System. The pump is removable so it can be used at main sink.

Mobile Cart Features

Optional
FlexiPump
™ Independent
Flushing System

Tools integrated
for pre-cleaning

Large inner well
that
accommodates
Endoscopes

PureStation
™ Sink Insert

Polypropylene (softer
material) reduces scope
damage.

Transport covers
available


Remove entire
insert and place
into existing sink



Versatile Pump



Flushing pump-
easily removes
so it can be
used at main
sinks.



FlexiPump can be
hung on the wall,
used on
countertop, or
hung from the sink
well'




Application for cart



- Short term displacement of processing
- Bring it out at peak processing times, put it in corner when not being used
- Sinks are too small for flexible scopes (coil scope to tight causing damage)
- Transportation- Lids available for transport

Draining Options

Ball valve drain



**Transfer
Pump**

