

PureSteel™ XI Reprocessing Sink

A small footprint pre-cleans larger devices.

The compact PureSteel XI Reprocessing Sink has a small footprint but allows a department to pre-clean a variety of instruments, including larger minimally invasive and robotic surgical devices and accessories. The sink is height-adjustable to bring pre-cleaning functions to an ergonomic height for the technician or nurse, and its basin depth is designed to eliminate deep bending. A 45-degree wrist rest is integrated on the front edge for comfortable, ergonomic positioning. This multifunctional unit can be mobile or stationary, and can be dedicated to specific devices or used to add pre-cleaning capacity to a small department. The XI sink is delivered with an accessory package to enhance ergonomics and productivity. Accessory package items may include: the FlexiPump™ Independent Flushing System, a stainless steel back wall to hold accessories, a custom-fit sink liner to protect instruments as they are processed, and one 2-gallon rinse and solution tank.



pure
PROCESSING

For more information on the PureSteel™ XI Reprocessing Sink, or to see a video demo:

VISIT www.pure-processing.com OR CALL 877.718.6868