


PureSteel™ Halogen Magnifier Task Light

A well-focused spot light for concentrated, visual inspection. This task light is water resistant for use in the decontamination room, and can be installed on desktops, walls, or back wall panels. Task lights are also ideal for prep and pack stations, or near sink stations to add light during pre-cleaning.



TECHNICAL SPECIFICATIONS

- Supply Voltage 230V/12V, 12V, 24V
- Light Source Halogen 12V/24V/20W
- Lifetime, light source 2,000 hours
- Arm Length 27 inches
- Color Black
- Color Temperature 3000K
- IP rating IP20
- Power Switch Placed in base
- Base Equipped with four mounting holes
- Cable Length 6.5 feet

pure
PROCESSING

TO SEE OUR VIDEO DEMONSTRATION GO TO
www.pure-processing.com/resources/videos/

For more information on the PureSteel
Halogen Magnifier Task Light, call:

877.718.6868

PureSteel™ LED Task Light

A high-intensity task light that uses less energy and enhances visual inspection capabilities. This task light is water resistant for use in the decontamination room, and can be installed on desktops, walls, or back wall panels. Task lights are also ideal for prep and pack stations, or near sink stations to add light during pre-cleaning.



27" flexible, vibration-resistant arm for precise control

LED intensity comparable to a 50 watt halogen light

Features on/off switch and dimmable lighting

Water resistance to work in moist environments

TECHNICAL SPECIFICATIONS

- Input 12-24V AC/DC; 50/60Hz (100-240VAC plugin AC/DC adaptor available, IP44)
- Power consumption 6W
- Arm length 27 inches
- Color Black
- IP rating IP65 (low voltage models)
- Beam angle 10° std, 25° optional
- Light output 270Lm, 36,000Lux at 50 cm, 10 degrees

pure
PROCESSING

TO SEE OUR VIDEO DEMONSTRATION GO TO
www.pure-processing.com/resources/videos/

For more information on the PureSteel LED Task Light, call:

877.718.6868