

Micro-Scientific® Goo Gone®

Topical and Surface Adhesive Remover



Micro-Scientific Goo Gone Topical and Surface Adhesive Remover removes tape and adhesive completely and safely from instruments, equipment and containers.

- Removes tape and adhesives quickly and effortlessly
- Does not harm skin while allowing removal of tape and adhesives without discomfort
- Product is safe and compatible with all metals

Product Features

Safe: Product is safe and compatible with all metals.

Rapid and effective: Removes tape and adhesives quickly and effortlessly.

Highly concentrated: Low cost per use.

Bacteriostatic: If store properly, not favorable to bacterial growth.

Environmentally-friendly: Biodegradable, phosphate-free formula.

Directions for Use

Saturate cotton swab or gauze with Micro-Scientific Goo Gone Topical and Surface Adhesive Remover and apply it to the edge of the tape or adhesive. Pull tape gently and keep applying Adhesive Remover. After using Adhesive Remover on instruments, equipment or containers, follow complete processing procedures again before proceeding with disinfection or sterilization.

pure
PROCESSING

To learn more, visit:
www.pure-processing.com

For more information, call:

877.718.6868