


Pure LandingPad™ Surface Liner

- ◆ Protects delicate devices and cushions surfaces where expensive endoscopes and other devices are staged for cleaning
- ◆ A ¾ inch thick nitrile rubber foam blend resists hospital-grade chemicals and disinfectants
- ◆ Won't absorb liquids
- ◆ Perforations facilitate the removal of liquids and debris from work surfaces
- ◆ For use on counter tops or staging areas
- ◆ 18X24 size



Pure LandingPad™ Silicone Sink Liner

- ◆ Place at the bottom of sink basins or counter tops to protect delicate instruments
- ◆ Non-toxic, eco-friendly, anti-static material
- ◆ Medical grade silicone
- ◆ 20x20 size
- ◆ Can be cut to size down from 20"x20" size



Protecting
Precious Instruments