



Pure LandingPad™ Surface & Silicone Liners – Care & Cleaning

Construction – Surface Pads

The Pure LandingPad™ Surface Liners are made of a PVC/nitrile rubber foam blend that does not absorb liquids.

Construction – Silicone Liners

The silicone liners are made of 100% medical grade silicone.

Cleaning

The Pure LandingPad™ Surface & Silicone Liners can be cleaned with a ready-to-use environmental disinfectant or in a cart washer on a setting that does not exceed 185F. Extended exposure to disinfectants with bleach may cause discoloration over time.

