

# PureSteel Healthcare Reprocessing Sinks

Pegboard Panels & Accessories



**Flexible**  
to meet your needs.

**Built**  
with instrument needs in mind.

**Inspired & designed**  
by our customers.

**Durable & cleanable**  
to meet demands.

**pure**  
PROCESSING

## Lumen Flushing

When decontaminating devices with lumens, personnel should... soak and flush the lumen according to the manufacturer's written IFU... Thoroughly flushing lumens helps ensure complete surface contact with the solution. If a brush is too large, it will not fit into the lumen; if it is too small, it will not have complete contact with the lumen walls and, consequently, will not clean them thoroughly.

– ANSI/AAMI ST79: 2017, 7.6.2



### FlexiPump Independent Flushing System

- Assures copious flushing of internal channels and lumens
- Reduces repetitive motion injuries from traditional flushing methods



## Basin Accessories



### PureStation™ Sink Inserts

- Raises working level in deep, pre-existing sink basins to improve employee ergonomics
- Features a temperature strip, volume markings, chemical mixing well, and overflow prevention handles
- Dedicated Eye Sink serves as a designated cleaning area for ocular instruments



NEW

### Staging Panel

- Create optional and flexible arrangements of additional counter space when staging room is required



NEW

### Splash Guard

- Provides a clear shield to cover a sink basin, reducing a user's contact with aerosols and contaminated reprocessing fluids
- Guards can be lifted when not in use



NEW

### Two Gallon Enzyme and Rinse Tanks

- Help separate reprocessing solutions from sink basins during the manual cleaning process
- Comes with quantity of 2, gallon containers and a bracket



NEW

### Basin Lights

- Integrated basin lights assist technicians with leak testing and visual inspection




### Silicone Surface Liner

- 20"x20" silicone basin liner protects delicate instruments
- For use in sink basins or counter tops



### Wrist Rests

- Removes pinch points and strain on users' wrists
- Diverts contaminated liquids away from technicians



### Volume Indicator Plate

- Removable, cleanable alternative to sink basin stickers
- Custom made to any sink basin size



NEW

### Workflow Plate

- Plates indicate soak, wash, rinse, and final rinse workflows
- Custom made to fit any sink basin



## Pegboard Accessories



### Pegboard Panels

- Individual panels for back wall organization that can be installed on any available wall



### Gallon Holder

- Keep gallon bottles organized on a back wall
- Holds 1-3 bottles



### Glove Box Holder

- Holder provides an area to place dispensable boxes of disposable gloves
- Custom sizes available



### Automated Endoscope Reprocessing (AER) Rack

- Plug adapters in close to the AER
- Can accommodate up to 11 tubes



NEW

### Bluetooth Speaker

- Connects by bluetooth to phones, computers
- Great for reprocessing sinks, tables



NEW

### Veriscan Leak Tester Shelf

- Shelf allows for storage of Veriscan Leak Tester off of countertop




### Small Leak Tester Shelf

- Shelf designed to support and mount smaller leak testers; great for Olympus-style leak testers



### Wire Baskets

- Organize reprocessing equipment and supplies onto back walls at ergonomic heights
- Multiple sizes available



### 32oz Bottle Holder

- Keeps bottles organized on a back wall
- Holds 1-3 bottles



### Waterproof Timer

- Submersible, waterproof timer for soaking stations and sinks
- Highly visible LED digits
- Includes bracket for installation



### Touchless Auto Dispenser

- Automated, hands-free enzymatic dispenser
- Provides dosage readout and temperature gauge



### Brush Bath and Shelf

- Allows instruments to soak on a back wall
- Built-in drain system



### Brush Holder & Organizer

- Organizes cleaning brushes at an easy to reach, ergonomic height
- Accommodates up to 8 brushes



### Pegboard Hooks

- Hang accessories from pegboards
- Hardware-free installation for rearrangeable layouts

## Faucets & Drains

NEW



### Rinse Sprayer

- Provides consistent and copious rinsing of instrument sets and trays using treated water

NEW



### Deionized Faucet

- Provides treated water
- Comes with a removable tip for attachment of accessories, such as a spray gun
- Handle is marked with a white DI designation per DI water applications



### Standard Faucet

- Wall mounted sink faucet
- Features 8" swing spout and lever handles



### Telescoping Drain Kit

- Pre-assembled drain systems for PureSteel sinks
- Plumbers only required for final connections and tightening



### Sprayer Faucet

- Add-on faucet features 44" flexible stainless steel hose, 10" swing nozzle

## Ergonomics

Ergonomic factors affecting worker safety and comfort should be considered when designing work spaces in the decontamination area, including...

c) anti-fatigue mats, constructed of materials capable of withstanding frequent cleaning, in areas where prolonged standing is required.

– ANSI/AAMI ST79:2017, 3.3.6.1.3



### GelPro Anti-Fatigue Mats

- Floor mats for sterile processing, operation, and surgical spaces



### GelPro Workflow Basin Mats

- Label sinks with appropriate reprocessing stages: soak, rinse, final rinse
- Cart washer safe

## Visual Inspection

Inspection using enhanced visualization tools such as lighted magnification and video borescopes might identify residues not observable by the unaided eye.

– ANSI/AAMI ST79: 2017, 7.6.4.5



### Magnification Task Light

- 1.8-3.0x magnification
- On/off switch and dimmable lighting
- 39" arm length



### Halogen Magnifier Task Light

- 2x magnifier
- Targeted illumination with 27" flexible goose-neck arm
- Water resistant



### Light Hood

- Durable light hood provides high-quality lighting at the sink
- Available in 4' and 2' lengths

pure  
PROCESSING

PHONE 877.718.6868 WEB pure-processing.com