



# PureSteel™ Halogen Magnifier Task Light

A well-focused spot light for concentrated, visual inspection. This task light is water resistant for use in the decontamination room, and can be installed on desktops, walls, or back wall panels. Task lights are also ideal for prep and pack stations, or near sink stations to add light during pre-cleaning.



## Lux rating - 0.5 m distance



## Technical Specifications

- ◆ Supply Voltage 230V/12V, 12V, 24V
- ◆ Light Source Halogen 12V/24V/20W
- ◆ Lifetime, light source 2,000 hours
- ◆ Arm Length 27 inches
- ◆ Color Black
- ◆ Color Temperature 3000K
- ◆ IP rating IP20
- ◆ Power Switch Placed in base
- ◆ Base Equipped with four mounting holes
- ◆ Cable Length 6.5 feet