



PureSteel™ Task Lights – Care and Cleaning

Manual Cleaning

Task lights should be disinfected routinely with intermediate level disinfectant sprays or wipes. Task lights should be turned off before being disinfected. Do not spray disinfectants directly onto bulbs. Follow the instructions given by your disinfectant manufacturer, as contact times may vary. It is recommended to use disinfectants low in bleach or alcohol when cleaning task lights. High concentrates of bleach or alcohol can cause fogginess or a dulling effect on magnifier lenses. A list of low bleach and low alcohol disinfectant options is provided below.

Product Name	Product ID
Spec Clean™ 24 oz Bottle	PSC240-1
Spec Clean™ Gallon Bottles	PSC128-1