

PureStation[®] Endoscopy Sink Insert

The PureStation Endoscopy Sink Insert's interior LED lighting helps technicians easily identify bubbles during leaking testing processes. The PureStation insert improves efficiency by providing convenient ways to monitor and mix proper enzymatic solution concentrations. Each basin includes an interior LED light, a replaceable thermometer strip, water volume indicators and a chemistry mixing well to make enzymatic use simpler and more consistent.

- Ergonomic solution for deep, utility style sinks
- The PureStation Endoscopy Sink Insert can be used in deep sinks, counters, and in mobile carts
- Lightweight for easy removal for when deep sink basins are required
- Cleaning tool caddy helps assure appropriate draining and drying of scrub brushes and device cleaning tools
- Chemical reservoir facilitates precise measuring of instrument cleaning chemistries

Specially designed catch drain system prevents small pieces from being lost down the drain

Easy-to-read thermometer band indicates the ideal temperature range for enzymatic detergents

Built-in handles facilitate carrying and provide an overflow prevention outlet



Interior LED light provides increased illumination and assists endoscopy departments identify bubbles during leak testing

Integrated watermarks on inner basin help achieve proper chemical-to-water

pure
PROCESSING

PHONE 877.718.6868 WEB www.pure-processing.com