

**POWERFUL · SIMPLE · EFFECTIVE**

# PureChannel™ Flush Instrument Channel Cleaner and Foam

PureChannel™ Flush is a triple-enzyme formula that is specifically designed to help remove the typical organic soils that stick to hard-to-clean internal channels and external surfaces after patient use, whether it's used as a diluted liquid cleaner or a concentrated foam application.



- Double-duty formula (liquid and foam) eliminates the need for two different products
  - Single dilution rate of 1/4 oz. per gallon helps eliminate confusion, inconsistency and waste
  - Concentrated foam breaks down soils, maintains moisture and facilitates cleaning
  - Regular use can help reduce the need for prolonged and intensive manual cleaning
  - Thoroughly tested and determined to be equally effective in a wide range of water conditions
  - Available in three ergonomic sizes: a one-gallon jug with handle, a lightweight 32-ounce bottle, and a 2oz foam bottle
- \* *Optional spray and foam nozzles sold separately.*

**pure**  
PROCESSING

For more information about our PureChannel™ Flush Instrument Channel Cleaner & Foam, please call: 877.718.6868 or, visit: [www.pure-processing.com](http://www.pure-processing.com)