



PureStation™ Mobile Cart

877.718.6868 | pure-processing.com



PureStation Mobile Cart



OSHA Central Supply,
“use carts with large, low rolling, low resistance wheels, that can easily roll over mixed flooring as well as gaps between elevators and hallways.”

The PureStation Mobile Cart creates another work station and is available with or without FlexiPump™ Independent Flushing System. The pump is removable so it can be used at main sink.

Mobile Cart Features



- Optional FlexiPump Independent Flushing System
- PureStation Sink Insert
- Remove entire insert and place into existing sink
- Tools integrated for pre-cleaning
- Polypropylene (softer material) reduces scope damage
- Large inner well that accommodates endoscopes
- Transport covers available

Versatile FlexiPump




- Flushing pump-easily removes so it can be used at main sinks.
- FlexiPump can be hung on the wall, used on countertop, or hung on a back wall panel.

Application for cart



- Short term displacement of processing.
- Bring it out at peak processing times; put it in corner when not being used.
- Sinks are too small for flexible scopes (coil scope to tight causing damage).
- Transportation-lids available for transport.

Draining options



Ball valve drain



Transfer Pump