


# PureSoak™ Timer

The PureSoak Timer allows departments to easily verify that timed reprocessing functions are being followed. Each timer can be programmed by either the seconds or minutes, and counts down until an alarm goes off, indicating the programmed time has elapsed. The large display makes it easy to see from across the room. Alarms are easily heard in even the busiest decontamination departments.



## Features

- ◆ Water resistant timer with large display to easily see across the department
- ◆ Programmable times that allow users to set times by either seconds or minutes
- ◆ Alarm feature when programmed times elapse
- ◆ Easy-to-use
- ◆ Great for reprocessing sinks or soaking areas
- ◆ Dimensions: 9.75"L x 4.75"H x 3"D

## Technical Specifications

- ◆ **Display:** Four digit bar; high intensity; seven segment; red LED display
- ◆ **Case:** All aluminum enclosure
- ◆ **Weight:** 3 lbs
- ◆ **Power Source:** 2 Watts; 90 - 240VAC 50/60 Hz, 6 ft 3-Prong power cord
- ◆ **Operating Temp:** 0° to 49° Celsius (32° to 120° Fahrenheit)
- ◆ **Case Finish:** Finger Print Resistant, Tough, Black Textured Powder Coating