

PureSteel™ Ergonomic WorkStation

PATENT PENDING

Meet instrument packaging and ergonomic standards

The PureSteel Ergonomic Workstation is designed to enhance the packaging of instruments. By providing illumination for the inspection of perforations in sterile wrap materials, staff can have more control over quality. Combined with an ergonomic workspace, worker safety is enhanced and productivity is increased. An array of accessories are also available to accommodate the specific needs of each department.

- Table top light provides illumination for inspection of sterile wrap materials
- Height-adjustability allows technicians to perform tasks while standing or seated
- Ergonomic-friendly drawer, shelving and accessories assist in diminishing repetitive task fatigue and injury
- Numerous accessories such as wire baskets, storage bins and electronic monitors are available
- Wide counter space provides more room for larger instruments and wraps
- Equipped with eight electrical outlets and underside cord management
- Casters allow for mobility
- Optional Organizational Rail organizes peel pouch rolls, stickers, and other sterilization materials easy dispensing and filling
- Optional Peel Pack Rack Shelf vertically stores and organizes pre-cut peel pouch packets on a back wall shelf



Meet compliance:

“When using woven wraps, sterile processing personnel should ... inspect wraps before each use for holes, worn spots, and stains.”

“When choosing and using non-woven wraps, sterile personnel should ...inspect wrap to ensure that it is free of defects that could have an adverse effect on the performance of the material.”

ANSI/AAMI ST79:2017, 9.1, Packaging Standards

Available in 48", 60", 72" and 78" lengths.



pure
PROCESSING

For more information:

VISIT www.pure-processing.com OR CALL 877.718.6868