

# PureSteel™

## Ergonomic Workstation and Wrap Inspection Tables




**pure**  
PROCESSING

PHONE 877.718.6868 WEB [www.pure-processing.com](http://www.pure-processing.com)

# PureSteel™ Workstation and Wrap Inspection Tables

Engineered by experts for understand the reprocessing space, PureSteel Ergonomic WorkStation and Wrap Inspection Tables are designed to enhance the inspection and packaging of instruments. By providing illumination for the inspection of instruments and perforations in sterile wrap materials, staff can have more control over quality. Height-adjustability and enhanced lighting help to improve staff safety and increase productivity. An array of accessories are also available to accommodate the specific needs of each department.

## PureSteel Ergonomic Workstation



To see a video demonstration, go to [www.pure-processing/resources/videos/](http://www.pure-processing/resources/videos/)

# PureSteel™ Packaging and Assembly Solutions

## PureSteel Wrap Inspection Table and Mini Wrap Inspection Table

The PureSteel Wrap Inspection Table is designed to organize and enhance the packaging and assembly of instruments. Combined with a back wall rack in which to hang wraps, the table top light provides illumination for the inspection of perforations in sterile wrap materials.

The PureSteel Mini Inspection Table is designed to provide sterile wrap inspection abilities in smaller spaces. The mobile table features an illuminated table top and lower shelf for storage of materials.



**pure**  
PROCESSING

PHONE 877.718.6868 WEB [www.pure-processing.com](http://www.pure-processing.com)