

# PureSteel™ Wrap Inspection Table

## Organize and inspect packaging and assembly wraps

The PureSteel Wrap Inspection Table is designed to organize and enhance the inspection packaging and assembly of instruments. Combined with a back wall rack in which to hang wraps, the table top light provides illumination for the inspection of instruments and perforations in sterile wrap materials.

- Table top light provides illumination for inspection of instruments and sterile wrap materials
- Back-wall rack allows for easy organization of wraps
- Dimmable switch controls lighting when inspecting wraps of different thickness
- Wide counter space provides room for larger instruments and wraps
- Bottom shelf allows for storage of materials
- Casters allow for mobility
- Available in 48", 60", 72" and 78" lengths. Custom sizes available



### Meet compliance:

"When using woven wraps, sterile processing personnel should... inspect wraps before each use for holes, worn spots, and stains."

- ANSI/AAMI ST79:2017, 9.1, Packaging Standards

# PureSteel™ Mini Inspection Table

The PureSteel Mini Inspection Table is designed to provide sterile wrap inspection abilities in smaller spaces. The mobile table features an illuminated table top and lower shelf for storage of materials.

- Table top light assists in visual inspection of instruments and wrap materials
- Small footprint
- Casters allow for mobility
- Bottom shelf creates space for storage of materials
- Dimmable switch controls lighting when inspecting wraps of different thickness



**pure**  
PROCESSING

For more information:

VISIT [www.pure-processing.com](http://www.pure-processing.com) OR CALL 877.718.6868