

# PureSteel™ Endoscopy Sink

The PureSteel™ Endoscopy Sink creates a dedicated space for endoscopy departments to soak, flush and rinse instruments per guidelines. The sink comes equipped with an LED-lit basin that assists technicians with leak testing. An integrated temperature gauge monitors enzymatic temperatures to meet IFU requirements.



- LED-lit basin facilitates clear visualization for leak testing
- Temperature gauge indicates enzymatic solution temperatures
- Solid back wall panel is easy to clean (pegboard available)
- Electric push button allows precise, personal height adjustment
- Ergonomic wrist edges remove pinch points
- Storage shelf below sink basins raises and lowers to allow mopping and cleaning, and easy access to stored items
- Integrated flushing systems automate the process for consistent, effective precleaning
- Meets all standards, guidelines and recommendations

*\*Standard and custom sizes available*

**pure**  
PROCESSING

For more information:

VISIT [www.pure-processing.com](http://www.pure-processing.com) OR CALL 877.718.6868