

# FlexiPump™ Independent Flushing System

Value based on increased processing capacity and optimized safety and compliance

---



All types of reusable lumened devices, including rigid and flexible scopes, robotic instruments, minimally invasive surgical devices, ocular devices and suction tips, must be flushed thoroughly and consistently to assure complete removal of debris and soils during pre-cleaning. Traditional manual flushing processes with syringes are subject to inconsistency and can be a cause of repetitive motion injury. The FlexiPump™ Independent Flushing System is a cost-effective automated alternative that helps assure compliant, high quality flushing, provides much greater workflow capacity, and supports worker and patient safety

## A huge capacity booster

Imagine a tech standing at the sink manually flushing a single duodenoscope channel over and over again with her 60 ml syringe. Now, imagine a department with three FlexiPump™ systems at that same sink (two on the wall and one on the counter), each pump with three tubes connected to a device, automatically flushing nine lumens at the same time, each in 60 seconds! It's not hard to see that having these compact, cost-effective systems at your sink can save time and greatly reduce manual labor, which can free up your technicians for other duties in the workflow.

## Improves reprocessing quality systems

National and international regulating bodies, such as the US Food and Drug Administration, the Joint Commission, AAMI, IAHCSSM, ECRI Institute and others, have been calling for more robust reprocessing quality systems in the wake of infection outbreaks around the world. They have flagged the need for standardized, consistent and repeatable processes, which can be enhanced by using automated systems where appropriate. The FlexiPump™ system injects consistent, repeatable flushing into the pre-cleaning process, with simple, intuitive operation that is easy to train and a system that's easy to maintain. It adds tremendous value for enhancing your department's reprocessing quality system.

## Non-disruptive installation

Any department can enjoy the benefits of using FlexiPump™ systems. It does not matter how small or how old your department is, you can easily find a spot for several FlexiPump™ systems near your pre-cleaning sink, and get going on making things better for your patients, your technicians, your reprocessing quality and your bottom line.



130 East Saint Charles Road, Unit C, Carol Stream, IL 60188 P 877.718.6868

[www.pure-processing.com](http://www.pure-processing.com)

FlexiPump™ is a trademark of Pure Processing LLC. © 2015, Pure Processing LLC. All rights reserved.

# FlexiPump™ Independent Flushing System

FlexiPump™ Features	Key Outcomes	Savings/Value to You
<b>Automated pump system</b>	<p>Significantly decreases manual labor and time required to adequately flush any lumened device</p> <p>Increases productivity and throughput</p> <p>Promotes worker and patient safety through a repeatable, consistent process</p>	<p>Reduced reprocessing cost</p> <p>Reduced technician time per device</p> <p>Opportunity to optimize personnel use in workflow</p> <p>Increased capacity and improved workflow</p> <p>Contributes to a quality system for flushing</p> <p>Fewer worker compensation claims and less lost work time</p>
<b>Versatile, non-disruptive installation</b>	<p>Requires no special installation or wiring, only a standard electrical outlet.</p> <p>Installs on a counter, on a wall, in upper edge of existing sink, or in/on Pure Station™ Medical Devices Pre-Cleaning System</p> <p>Simple nest and stick-on rail accessories for counter, wall and sink installation</p> <p>Can also be used with an optional Pure Station™ Medical Devices Pre-Cleaning System on a rolling cart</p>	<p>Does not add to department footprint</p> <p>Does not produce any remodeling cost</p> <p>Technicians can easily install the system without additional assistance</p> <p>Low-cost non-disruptive accessories customize maximum functionality for your facility</p>
<b>Flexible connectivity</b>	<p>System connects to any channeled device, with or without luer locks, with or without ports, and even devices without either luer locks or ports</p>	<p>System assures repeatable, consistent, quality flushing for your entire channeled device inventory</p>
<b>Easy operation</b>	<p>Four buttons, easy-to-learn functions</p> <p>Basic, 5-step maintenance to decontaminate the tube set and 4-step process to replace it</p>	<p>Plug-and-play design requires minimal training</p> <p>System is easy to maintain in optimal working order, which enhances the department's quality system</p>

*\* The FlexiPump™ Independent Flushing system is intended for pre-cleaning and flushing of lumened devices. It is not for final treatment and does not replace manual brushing and rinsing. Refer to each medical device manufacturer's instructions for use for complete pre-cleaning directions.*

