

# PureSteel™ XL Transport Carts

Easily transport and soak robotic instruments in InstruSafe® XL Transport Containers

Transporting robotic and larger surgical instruments in soaking containers can pose ergonomic and biohazard safety risks. The transport containers and carts provide a secure, dedicated space to horizontally soak instruments while in transit from the operating room to the decontamination area. The carts also serve as a sealed soaking station for facilities lacking adequate soaking space.

PureSteel™ XL Transport Container Cart

PureSteel™ XL Transport Container Tiered Cart

Organizational back wall allows for installation of the FlexiPump™ Independent Flushing System, PureSteel Waterproof Soak Timer and other accessories

Cart shelf is designed to facilitate container drains for safe, easy disposal

Shelving provides vertical storage for InstruSafe XL Transport Containers and container lids

Casters provide mobility

Provides a staging area for soaking up to four levels of soiled devices at a time

Each cart accommodates vertical storage for up to four InstruSafe XL Transport Containers



For more information:

VISIT [www.pure-processing.com](http://www.pure-processing.com) OR CALL 877.718.6868