

PureSteel™ Wrap Inspection Table

Organize & inspect packaging & assembly wraps

The PureSteel Wrap Inspection Table is designed to organize and enhance the packaging and assembly of instruments. Combined with a back wall rack in which to hang wraps, the table top light provides illumination for the inspection of instruments and perforations in sterile wrap materials.

- Table top light provides illumination for inspection of instruments and sterile wrap materials
- Backwall rack allows for easy organization of wraps
- Dimmable switch controls lighting when inspecting wraps of different thickness
- Wide counter space provides room for larger trays and wraps
- Bottom shelf allows for storage of materials
- Available in 60", 72" and 78" lengths. Custom sizes available



Ensure compliance:

"When using woven wraps, sterile processing personnel should... inspect wraps before each use for holes, worn spots, and stains."

- ANSI/AAMI ST79:2017, 9.1, Packaging Standards



PureSteel™ Mini Inspection Table

The PureSteel Mini Inspection Table is designed to provide sterile wrap inspection abilities in smaller spaces. The mobile table features an illuminated table top and lower shelf for storage of materials.

- Table top light assists in visual inspection of instruments and wrap materials
- Small footprint
- Bottom shelf creates space for storage of materials
- Dimmable switch controls lighting when inspecting wraps of different thickness

PureSteel™ Mobile Scale Cart

The PureSteel™ Mobile Scale Cart offers a secure space for weighing trays with mobility features like casters and a rechargeable battery. Its built-in display allows users to handle, record, and print tray weights directly from the scale for convenience.

- The built-in scale lets users easily weigh trays and instruments
- Under-shelf allows for additional storage



PureSteel™ Ergonomic WorkStation and Wrap Inspection Tables



pure
PROCESSING

PHONE 877.718.6868 WEB www.pure-processing.com

pure
PROCESSING

PHONE 877.718.6868 WEB www.pure-processing.com

PureSteel™ Ergonomic WorkStations & Wrap Inspection Tables

Engineered by experts who understand the reprocessing space, PureSteel Ergonomic WorkStations and Wrap Inspection Tables are designed to enhance the inspection and packaging of instruments. By providing illumination for the inspection of instruments and perforations in sterile wrap materials, staff can have more control over quality. Height-adjustability, a backwall pegboard and customizable built-in features help to improve staff safety and increase productivity. An array of accessories are also available to accommodate the specific needs of each department.

- Optional tabletop light provides illumination for inspection of sterile wrap materials and surgical instruments
- Height-adjustability allows technicians to perform tasks while standing or seated
- Wide catalog of accessories, from monitors to shelving, to tailor tables to their needs and use
- Optional, fully customizable back walls organize accessories, power cords and computers at the most ergonomic height for technicians; adjust up & down as the table height adjusts
- Wide counter space provides more room for larger instruments and wraps
- Equipped with eight electrical outlets and cord management
- Standard sizes: 60", 72" and 78" lengths
- Custom sizes available



To see a video demonstration, go to www.pure-processing/resources/videos/