




PureSteel™ Magnifying LED Task Light

Inspect instruments with ease



The PureSteel Magnifying LED Task Light makes inspection of instruments and devices easier with 3.85x magnifying power and selective LED lighting. PureSteel Magnifying LED Task Lights can be installed on tabletops and back wall panels. PureSteel LED Task Lights are also ideal for meeting packaging and assembly guidelines and standards.



Errors and technical changes excepted. Images can vary from original.
Photometric and electrical data may vary due to component tolerances.

Technical Features

- ◆ Outfitted with 48 LED lights
- ◆ Multiple LED sections can be switched on or off and are dimmable
- ◆ Exchangeable magnification lenses easily remove for cleaning
- ◆ Features 39" stable arm, sealed joints, robust aluminum and high-performance plastic
- ◆ Voltage 100-240V, 50/60Hz
- ◆ Operation touch-panel on head
- ◆ Magnification 3.5 diopters or 3.5d + 8d
- ◆ Color Temperature 4000K
- ◆ Color Rendering Index (CRI) 90
- ◆ Luminaire Material anodized aluminum
- ◆ Lens Material acrylic
- ◆ Power Consumption 14W
- ◆ Protection Rating IP20
- ◆ Certification CE, ETL/cETL
- ◆ Cord/Plug 12' cord w/ US plug