

# Micro-Cide™ 28 HLD

## 3% Glutaraldehyde Reusable Sterilizing and High-Level Disinfection Solution



Micro-Scientific® Micro-Cide™ 28 HLD is a 3% glutaraldehyde solution with attached powder activator. Once activated, Micro-Scientific® Micro-Cide™ 28 HLD can be used for 28 days to provide high-level disinfection or sterilization, providing the solution is above the 1.8% minimum effective concentration and is used according to directions for use.

- Provides a way to effectively sterilize or achieve high-level disinfection of heat sensitive reusable medical devices
- Achieves high-level disinfection of instruments in 25 minutes, providing directions for use are followed
- Sterilizes instruments in 10 hours, providing directions for use are followed



To learn more, visit: [www.pure-processing.com](http://www.pure-processing.com)

For more information, call: 877.718.6868