


Prep and Pack Table Comparison Guide

No two prep and pack tables are created equal. From time and motion for worker productivity, to ergonomics and its impact on staff, to organization and its impact on turnover time and compliance, assembly tables require careful planning and consideration to ensure they meet each department's individual needs and objectives.



Pegboard and wrap rung options organize both traditional table accessories, and integrate wrap storage at the point of use.



EZ Slide Scale weighs surgical instrument sets and trays, compliant with AAMI and AORN.



Optional, integrated table lights enhance both wrap and instrument inspection.

FEATURE	BENEFIT & IMPORTANCE	PURE PROCESSING	COMPETITOR
Stainless grade	Durability and longevity of equipment	Made of 304 stainless steel for all day, heavy instrument reprocessing	Variety of grades & quality
Custom options	Work within constraints you can't change, or fully optimize existing or new layouts & floorplans	Depth, height, length all customizable	Limited options only
Installation	Arrives pre-installed with accessories to reduce install time	Comes standard, all accessories pre-installed before shipment for easy integration	Not available
Height-adjustable lifters	Increased ergonomics for all staff; OSHA, AORN, ANSI/AAMI ST79 compliant	3 year warranty on electrical lifters; reducing downtime and increasing ease-of-serviceability	Variety of lifters
Tabletop light	Technicians can inspect wrap for quality before use; no need to hold wrap up to dim ceiling lights Great for inspecting small instruments for functional integrity	Optional	Warranty & servicing requirements unknown

FEATURE	BENEFIT & IMPORTANCE	PURE PROCESSING	COMPETITOR
Height-adjustable shelf	Reduces back bending for equipment stored under the table; allows total visibility under the table and ease of cleaning	Available, optional	Not available
Electrical	More power for growing equipment needs on tables; doesn't require facilities to add outlets	8 integrated GFCI electrical outlets for future capacity; 1 outlet required per table	Power strips typical
Power management	Wires are safely organized behind and under the tables, away from heavy trays and counters	Standard; part of the Installation process	Done by hospital facilities staff, typically
Pull-out counter	Countertop which can be pulled out when additional space is needed	Available, optional	Not available
EZ Slide Scale	Slide scale allows you to weigh trays when needed	Available, optional	Not available
Undertable shelves	Multiple, suspended shelf options dependent on the need; ideal for storing printers or PC towers Adjust with height-adjustable tables, and keep accessories at more ergonomic heights	Available, multiple options	Limited options only
Heat Sealer Shelf	Mount heavy & bulky heat sealers on the underside of assembly tables. Shelf is ergonomic and collapsible, providing additional space.	Available, optional	Not available
Computer, IT installations	Upgrade IT systems, or use existing computers with enhanced mounting arms	Full systems available; multiple mounting options available	Not available
Footrest	Add ergonomic comfort for staff	Available, optional	Limited options only
Light hoods	LED light hoods provide long-term lighting to increase visibility at the table	Available, optional	Available
Etched table ruler	Ruler etched into table surfaces to assist with inspection and assembly	Available, optional in multiple lengths and configurations	Not available
Future additions	Future accessories & developments easily integrate into existing tables; table designs change as instrument reprocessing evolves	Regularly updated with new accessories and features	Not available

