

PureSteel™ Wrap Inspection Table

Organize and inspect packaging and assembly wraps

The PureSteel Wrap Inspection Table is designed to organize and enhance the packaging and assembly of instruments. Combined with a back wall rack in which to hang wraps, the table top light provides illumination for the inspection of perforations in sterile wrap materials.

- Table top light provides illumination for inspection of sterile wrap materials
- Back-wall rack allows for easy organization of wraps
- Dimmable switch controls lighting when inspecting wraps of different thickness
- Wide counter space provides room for larger instruments and wraps
- Bottom shelf allows for storage of materials
- Casters allow for mobility
- Available in 60" length. Custom sizes available



Meet compliance:


"When using woven wraps, sterile processing personnel should... inspect wraps before each use for holes, worn spots, and stains."

- ANSI/AAMI ST79:2017, 9.1, Packaging Standards

PureSteel™ Mini Inspection Table

The PureSteel Mini Inspection Table is designed to provide sterile wrap inspection abilities in smaller spaces. The mobile table features an illuminated table top and lower shelf for storage of materials.

- Table top light assists in visual inspection of wrap materials
- Small footprint
- Casters allow for mobility
- Bottom shelf creates space for storage of materials
- Dimmable switch controls lighting when inspecting wraps of different thickness



pure
PROCESSING

For more information:

VISIT www.pure-processing.com OR CALL 877.718.6868