

PureSteel™ Signature Series Ergonomic WorkStations

The PureSteel™ Signature Series Ergonomic WorkStations are designed to enhance the packaging and inspection of surgical trays & instruments.

A fully-optimized ergonomic prep & pack table, the Signature Series workstations improve worker safety and productivity through simple, thoughtful innovations:

- ◆ A height-adjustable table & pegboard that keeps activities and tasks at the most ergonomic height, per OSHA recommendations.
- ◆ Overhead lighting meets ST79 luminosity requirements and makes detailed inspection easier.
- ◆ Pegboards organize a wide range of accessories and help wire power away from table surfaces; monitors, assembly tools, and accessories can be freed from the countertop to create more space for assembling large trays and packs.
- ◆ Electro-mechanical lifters, which require little to no maintenance or care in their lifetime.
- ◆ Highest quality, 304 stainless steel for durability & longevity.
- ◆ Arrives pre-installed with accessories and features; no downtime waiting for installation.
- ◆ Future pegboard accessories can easily integrate into existing tables; tables grow with the needs of your department.



Technical Specifications

Standard options only: 72"x32" & 60"x32"

Height-adjustable

The Signature Series Ergonomic WorkStations also come with the full line of PureSteel accessories, including:

- ◆ Optional built-in EZ Slide Drawer for storage.
- ◆ Optional built-in EZ Slide Scale to weight trays as they're being assembled.
- ◆ Optional EZ Slide Countertop for added assembly space.
- ◆ Task lights for detailed inspection (options up to 4x magnification available).
- ◆ Lighted shelving, peel pouch organizers, bins, baskets, and more, to organize tools and items.