

Optimize prep-and-pack in your space

The PureSteel Wrap Inspection Table helps organize and enhance prep and assembly functions. The back wall rack organizes wraps to speed retrieval during packing. The tabletop light facilitates instrument inspection and helps to identify perforations in sterile wraps.

The PureSteel Mini Inspection Table provides sterile wrap inspection abilities in smaller spaces. The mobile table features an illuminated tabletop and lower shelf for storing materials.



pure
PROCESSING


130 East St. Charles Road, Unit C
Carol Stream, IL 60188

For more information, please call or visit:

877.718.6868

www.pure-processing.com

Visit us on our social media sites!



PureSteel™ Ergonomic WorkStation

Designed to enhance instrument inspection and packaging

"The tabletop light "solves issues. Holes in the filter or drape is a problem, so the light is a time-saver and a frustration-saver. We can have faith that the drape will be in good shape. When Central Sterile (CS) inspects instruments before [assembling] trays, they need light to find bioburden. The light on top of the prep and pack table gives them extra light."

- Sue, Sterile Processing Department Manager

pure
PROCESSING

Workstations should optimize ergonomics and efficiency while also meeting standards

"The clean work area/room should include space for packaging materials and preparation supplies; computers, if used, along with computer accessories; magnifying lights; and processing tables, which should be made of nonporous materials (stainless steel), ergonomic, and preferably, height-adjustable"

- ANSI/AAMIST79:2017 3.3.6.2 Clean work area/room

- + **Height adjustability:** push-button, table height adjuster meets ergonomic guidelines.
- + **Height adjustable back wall** organizes equipment, accessories and tools at ergonomic heights and frees up tabletop work space.
- + **Accessories** such as magnified LED task lights, wire baskets, storage bins and electronic monitors help optimize workflows.
- + **Cord management** features eight electrical outlets and underside cord organization.
- + **Organizing rail and peel pack shelf** locates peel pouch rolls, stickers and other materials for easy dispensing.
- + Magnifying task lights provide up to 4x magnification
- + **Pull-out counter, drawer and PC shelf** options allow flexible workflow enhancements.
- + **Sizes:** 48", 60", 70", 72" and 78"
Custom sizes available.




PHONE 877.718.6868

Inspections are critical to assure patient safety

"When using woven wraps, sterile processing personnel should ...inspect wraps before each use for holes, worn spots, and stains. ...When choosing and using non-woven wraps, sterile personnel should ...inspect wrap to ensure that it is free of defects that could have an adverse effect on the performance of the material."

ANSI/AAMI ST79:2017, 9.1, Packaging Standards



- + Tabletop light provides dimmable illumination for inspecting instruments and sterile wraps.
- + Task lights and hood lights meet American National Standard A11.1-1965, R1970, Practice for Industrial Lighting to reduce eye strain and assist in aiding the accomplishment of detail-intense tasks.

WEB pure-processing.com