

# PureSteel™ Mobile Scale Cart

The PureSteel™ Mobile Scale Cart offers a dedicated, secure space for weighing trays during assembly. Designed for mobility, this cart features casters and is equipped with a rechargeable battery, enabling multiple technicians to utilize it throughout the department without the need to be plugged into the wall. With the attached display, users can conveniently hold, and add tray weights directly from the scale.

*ANSI/AAMI, 4.3.5 Maximum Weight: The combined weight of the containment device, the instruments, and any accessories or wrappers shall not exceed 25 pounds when the containment device load is configured according to the manufacturer's written IFU.*



The built-in scale lets users easily weigh trays and instruments

Under-shelf allows for additional storage

Casters provide mobility

Tare button allows you to see the current change in weight when instruments are added or removed from a tray

Holds the weight reading of a tray without having tray present after initial read

Total button allows you to accumulate weight of multiple trays or instruments and gives a total



**pure**  
PROCESSING

PHONE 877.718.6868  
WEB [www.pure-processing.com](http://www.pure-processing.com)

