

PureSteel™ Signature Series Ergonomic WorkStations

The PureSteel™ Signature Series Ergonomic WorkStations are designed to enhance the packaging and inspection of surgical trays & instruments.

A fully-optimized ergonomic prep & pack table, the Signature Series workstations improve worker safety and productivity through simple, thoughtful innovations:

- A height-adjustable table & pegboard that keeps activities and tasks at the most ergonomic height, per OSHA recommendations.
- Overhead lighting meets ST79 luminosity requirements and makes detailed inspection easier.
- Pegboards organize a wide range of accessories and help wire power away from table surfaces; monitors, assembly tools, and accessories can be freed from the countertop to create more space for assembling large trays and packs.
- Electro-mechanical lifters, which require little to no maintenance or care in their lifetime.
- Highest quality, 304 stainless steel for durability & longevity.
- Arrives pre-installed with accessories and features; no downtime waiting for installation.
- Future pegboard accessories can easily integrate into existing tables; tables grow with the needs of your department.



Technical specifications

Standard options only: 72"x32" & 60"x32"
Height-adjustable



The Signature Series Ergonomic WorkStations also come with the full line of PureSteel accessories, including:

- Optional built-in EZ Slide Drawer for storage.
- Optional built-in EZ Slide Scale to weight trays as they're being assembled.
- Optional EZ Slide Countertop for added assembly space.
- Task lights for detailed inspection (options up to 4x magnification available).
- Lighted shelving, peel pouch organizers, bins, baskets, and more, to organize tools and items.