


APPLICATION

In the 1st Stage, the Robust Leading Brush dislodges debris; in the 2nd Stage, the Squeegee multi-disk wiper cleans the walls; and in the 3rd Stage, the Fine Bristled Tail Brush removes remaining debris and biofilm, ensuring a pristine finish before the final detergent flush.

DESCRIPTION

The SuperDART® ONE PASS Brush is a validated, three-stage, single-use endoscope cleaning device that efficiently eliminates debris and biofilm from channels, providing high-level cleaning results and time-saving advantages with reduced fatigue for nurses and technicians. Calibrated and color coded SmartSTOP™ feature reduces mistakes by preventing channel jams and detecting channel flaws.



The Selections Checked Below Apply to This Equipment

SuperDART® Part Number	Kit Part Number (includes a shorter, cylinder channel cleaning brush)	Catheter Length (cm)	Catheter Diameter (cm)	Units/box
CHANNEL SIZE RANGE 2.0 - 2.8MM Ideal for Cystoscopes, Bronchoscopes, Intubation Scopes...				
9293 <input type="checkbox"/>	9399 <input type="checkbox"/>	240	1.6	25
9253 <input type="checkbox"/>	9355 <input type="checkbox"/>	122	1.6	25
9233 <input type="checkbox"/>	9333 <input type="checkbox"/>	85	1.6	25
CHANNEL SIZE RANGE 2.8 - 4.2MM Ideal for Colonoscopes, Gastrosopes, Duodenoscopes...				
9292 <input type="checkbox"/>	9188 <input type="checkbox"/>	250	1.7	25

Facility

Contact Name

Facility Address

Representative Name

Contact E-mail

Representative E-mail

Date