

PureSteel™ Reprocessing Solutions



PureSteel™ Xi Healthcare Reprocessing Sink

- Height-adjustable for ergonomics
- Slanted wrist edges to prevent wrist strain & injury
- Small footprint
- Integrated with additional electrical outlets for capacity
- Creates a dedicated reprocessing area for instruments with unique IFU or requirements
- FlexiPump™ Independent Flushing System, Touchless Dose Pump, protective sink liners, and Waterproof Timers available

PureSteel™ Mobile Soaking Station

- Creates extra sink spaces for reprocessing
- Creates dedicated spaces for soaking, and is ideal for 30 minute robotic soaks per IFUs
- Accommodates robotic instruments & orthopedic sets
- Small, mobile footprint conserves space
- Great for flexible and rigid endoscopes



PureSteel™ Healthcare Reprocessing Sinks



PureSteel™ Healthcare Reprocessing Sinks

Engineered by sink experts who understand the reprocessing space, PureSteel™ Healthcare Reprocessing Sinks are the comprehensive stainless steel sink choice for reprocessing departments.




Newest Features



Paladin Evolution Rail System

- Easily rearrange accessories for flexible workflows.
- Customize shelf placement anywhere along the rail.



Copper Clean

Copper Clean® Anti-Microbial Wrist Rests provide continuous anti-microbial, self-sanitization properties in high-touch environments



PureSteel™ Staging Panel

Each Staging Panel can be placed onto a sink basin to create extended work surfaces for reprocessing.



PureSteel™ EZ Hinge™ Pegboard

- Optional, hinge pegboards help technicians and nurses clean behind pegboards or back walls
- In tight alcoves or installations, hinge pegboards provide control to move accessories, access hardware, utilities and more
- Requires 18"H backsplash



PureSteel™ UV LED Blacklight Light Hood

- UV blacklights to aid in hands-free detection of residual cement from orthopedic devices
- Can be used independently of regular overhead LED lighting