




2026 Product Catalog

◆ Raise the Standard



The Standard Has Been Raised - Our Look, Too

Reprocessing is the hidden backbone of safe healthcare. Every day, sterile processing and GI teams stand between patients and infection risk. Their work is essential—but often unseen, physically demanding, and under-supported. Pure Processing exists to champion you & your dedicated teams by developing and evolving solutions that protect you, your patients, and improve the operational and financial future of your facilities.

Why We Rebranded

Our new look wasn't created just to be different. At its core, our new look expresses our belief that great reprocessing elevates people, processes, and outcomes—and that your work deserves tools designed with that same standard of excellence.

The Story Behind our New Logo

Clarity + Excellence

The diamond accent symbolizes that moment when an instrument, workflow, or process is complete and ready for safe patient care. Our new diamond also represents our company's north star: you! Our customers are the guiding, focal point in everything we do.



Precision + Partnership

A bold "P" represents precision and partnership. It's a reminder that when we work together with our customers and take care of every detail with consistency and care, we achieve great things.

Our Promise to You

Every product we make and every process we improve exists for one reason: **to raise the standard of safety, efficiency and dignity in reprocessing.**

We are grateful for all your faith and dedication in our Mission these past 15 years. We look forward to being your Partnership for the next 15, and many more!

Table of Contents

PAGES
4-15



Decontamination Solutions

Effective, simple and validated technologies for standardized & high-quality results

PAGES
16-25



Healthcare Reprocessing Sinks & Accessories

Improving ergonomics, productivity and organization with custom solutions

PAGES
26-28



Mobile Solutions & Carts

Improving throughput and capacity in even the smallest departments


PAGES
29-33



Prep and Pack Tables & Accessories

Increasing productivity, organization and space with custom, detailed design

PAGES
34-38



Quality Assurance Technologies

Reducing repair costs, reprocessing errors, and creating a culture of high-quality outcomes



Decontamination Solutions



GSA Advantage!

Now on FSS contract! Find us on the GSA Advantage website.



FlexiPump™ Independent Flushing System

- Flush multiple instruments at once
- Adaptable to virtually all lumened instrumentation, including flexible and rigid scopes, robotics, laparoscopic, ocular, and more
- Provides hands-free, ergonomic manual cleaning
- Pressure relief valves help protect delicate instruments during flushing


QUANTITY	PRODUCT ID#
1	PSFPS-01

“


The FlexiPump made decontamination for cannulated instruments a lot easier for not just robotics, but for everything. That's why I'm not having to flush with syringes anymore. It made things a lot easier for staff and decontamination.

Matt, Central Sterile Manager

”



Syringe: *thoroughly* inconsistent, inefficient and ineffective



FlexiPump™ Independent Flushing System: *thoroughly* consistent, standardized and compliant



FlexiPump Brackets

- Comfortably nests one, two or three FlexiPump Systems onto a back wall or pegboard
- Made from stainless steel

QUANTITY	PRODUCT ID#
1 (single bracket)	FP-SNBR
1 (double bracket)	FP-DBBR


FlexiPump Irrigation Tube Sets

- Disposable tube sets for FlexiPump Independent Flushing System
- Available in 15 or 30 psi settings per device IFUs
- Available in 3 or 5 lead options
- Pre-assembled tube sets come with a new cap, rotor and ceramic weight and strainer with each tube set for convenience & easy installation
- Must be replaced every 30 days

Tube sets
should be
replaced every
30 days



PSI	QUANTITY	PRODUCT ID#
5 LEAD		
15 PSI	1 box of 6	FXITS-0515psi6-PH
15 PSI	6 pre-assembled	FXITS-0515psi6-PH-PREA
<i>For instruments smaller than 1mm in diameter only.</i>		



3 LEAD		
15 PSI	1 box of 6	FXITS-0315psi6-PH
15 PSI	6 pre-assembled)	FXITS-0315psi6-PH-PREA
30 PSI	1 box of 6	FXITS-0330psi6-PH
30 PSI	6 pre-assembled)	FXITS-0330psi6-PH-PREA



FlexiPump Accessory Kit

- Includes 50 XL Slip Tips, 25 End Caps, 20 Female to Female Luers, and 5 Cup Adapters
- Attaches to instruments not adaptable with a male luer lock on a syringe, including robotics, ocular instruments, and reamers

QUANTITY	PRODUCT ID#
1	ACC-KIT



FlexiPump Cap

QUANTITY

PRODUCT ID#

1 pack of 10

FLX-CAP-PH



FlexiPump Rotor

1 pack of 10

FLX-ROTOR-PH



End Cap

1 pack of 100

ENDCAPS-100



Rubber Cup Adapter

1 pack of 100

RCA-100



Female-to-Female

1 pack of 100

FFLL-100



XL Slip Tip

1 pack of 100

XLSLIP-100



TimeStrip® Time Indicators

- TimeStrip® Time Indicators offer a helpful visual reminder of service and maintenance schedules for consumable components

100 (1 month)

TM-1M

100 (3 month)

TM-3M



FlexiPump Ceramic Weight and Strainer

- Anchors the drop-down tube of FlexiPump tube set beneath water

1 pack of 10

FPCW-10



FlexiPump Power Supply

1

FPPS-01



PureChannel Flush™ Touch-less Auto Dispenser

- Automated, hands-free enzymatic dispenser to reduce chemical handling and prevent under- and over-dosing, while providing a dosage readout and temperature gauge



QUANTITY	PRODUCT ID#
1	PCFTAD-01
1 (replacement kit)	PCFTAD-RKIT



PureChannel Flush Instrument Channel Cleaner and Foam

- A triple enzyme surfactant formulation that is designed for soaking and flushing medical devices with hard-to-clean narrow channels
- Dilutes at 1/4 oz. per gallon
- Highly effective as both a liquid and foam
- Easy to use with the FlexiPump Independent Flushing System, in sonic cleaners, or for manual cleaning in a sink



1 case (12, 32 oz. bottles)	PCF-32oz
1 case (4, 1 gallon bottles)	PCF-08



Weiman Stainless Steel Cleaner

- Weiman Stainless Steel Cleaner shines stainless steel equipment with a streak-free shine, while repelling fingerprints and dust.

1 case (6, 22 oz. trigger bottles)	WSST
1 box (4 packs of 30 wipes)	WSSW
1 case (6, 17 oz. aerosol bottles)	WSSA



Spec Clean™ Ready-to-Use Disinfectant/Cleaner

- One step ready-to-use disinfectant/cleaner
- Safe for use on electronic panels, PureStation Sink Inserts, and FlexiPump tubing
- Sample kits per request*



1 case (15, 24 oz. bottles)	PSC240-1
1 case (4, 1 gallon bottles)	PSC128-1

NEW

SuperDART®
ONE PASS™ Validated

1st
STAGE


ROBUST LEADING BRUSH

2nd
STAGE

SQUEEGEE MULTI DISK WIPER

3rd
STAGE

FINE BRISTLED TAIL BRUSH



SuperDART® ONE PASS™ Validated w/SmartSTOP™

Single-Use Endoscope Channel Cleaning Device

- Clean endoscope channels in just one pass
- Easy-pull-through action with a robust leading brush, and second squeegee multi-disk wiper, and third fine bristled tail brush
- Three stage brush which completes the 3-pass cleaning process in just ONE PASS™
- SmartSTOP™ tip prevents accidental jamming of brush inside channels, reducing errors & repairs costs



LENGTH	QUANTITY	PRODUCT ID#
Channel Range: 2.0mm-2.8mm (Brush Only)		
240cm	25 - Box	9293
122cm	25 - Box	9253
85cm	25 - Box	9233
Channel Range: 2.0mm-2.8mm (Brush + Suction Cylinder Brush)		
240cm	25 - Box	9399
122cm	25 - Box	9355
85cm	25 - Box	9333
Channel Range: 2.8mm-4.2mm (Brush Only)		
250cm	25 - Box	9292
Channel Range: 2.8mm-4.2mm (Brush + Suction Cylinder Brush)		
250cm	25 - Box	9188

Channel Range: 2.0mm-2.8mm (Brush Only)

240cm	25 - Box	9293
122cm	25 - Box	9253
85cm	25 - Box	9233

Channel Range: 2.0mm-2.8mm (Brush + Suction Cylinder Brush)


240cm	25 - Box	9399
122cm	25 - Box	9355
85cm	25 - Box	9333

Channel Range: 2.8mm-4.2mm (Brush Only)

250cm	25 - Box	9292
-------	----------	------

Channel Range: 2.8mm-4.2mm (Brush + Suction Cylinder Brush)

250cm	25 - Box	9188
-------	----------	------





PureStation™ Sink Insert

GSA Advantage!

Now on FSS contract! Find us on the GSA Advantage website.

- Raises working level in deep, pre-existing sink basins to improve employee ergonomics.
- Includes integrated accessories for easier enzymatic maintenance: overflow prevention handles, volume markings, thermometer strip, a chemical dilution well, solution and brush caddies, and a drain filter



“Ergonomic factors affecting worker safety and comfort should be considered when designing work spaces in the decontamination area...”

ANSI/AAMI ST79:2017 3.3.6.1.3

”

QUANTITY	PRODUCT ID#
1 (16"x16")	PS1616
1 (16"x26")	PS1626



Small Device and Solution Caddies

- 2 solid solution containers and 1 small device container with perforations for easy draining

1 kit of 3 DSC-KIT



PureStation Thermometer Strip

1 pack of 10 THERMO-10



Drain Plug with Chain

- Tethers drain plug to sink insert

1 pack of 5 DPWCHAIN



PureStation™ Dedicated Eye Sink

- Serves as a separate designated cleaning area for channeled ocular surgical devices
- Includes special labels indicating insert for use only with ocular devices; reminder to keep ocular reprocessing separated from other specialties

PURE
PROCESSING™

pure-processing.com

EYE Instruments Only

Discard H₂O used to clean or rinse after each use. Follow device manufacturer's instructions for use.

(877) 718-6868




“Intraocular surgical instruments should be processed separately from general surgical instruments and equipment to reduce the potential for cross-contamination by material or residue from general surgical instruments.

ANSI/AAMI ST79:2017 Annex M, M.2.2.2

”



QUANTITY	PRODUCT ID#
1 (16"x16")	PS1616-EYE
1 (16"x26")	PS1626-EYE



Digital Thermometer

- Displays cleaning solution temperature
- Programmed to alarm when temperature is out of range

1

PS-DTHERMO



GelPro™ NewLife Eco-Pro Anti-Fatigue Mats

- Durable floor mats with a high-traction bottom, low-profile edges that prevent tripping, and easy-to-clean water-resistant matte surfaces
- Includes yellow stripe to indicate beveled edge

QUANTITY	PRODUCT ID#
1 (20" x 32")	GPAFM20x32
1 (20" x 48")	GPAFM20x48
1 (20" x 72")	GPAFM20x72
1 (24" x 36")	GPAFM24x36
1 (36" x 60")	GPAFM36x60



GelPro Workflow Mats

- Label sink basins with their appropriate reprocessing step, and keep tasks separated according to guideline recommendations
- Provide ergonomic relief from long periods of standing at the sink

1 (20" x 32"), Soak Mat	BGPAFM20x32-SK
1 (20" x 32"), Wash Mat	BGPAFM20x32-WS
1 (20" x 32"), Rinse Mat	BGPAFM20x32-RS
1 (20" x 32"), Final Rinse Mat	BGPAFM20x32-FS



Black Only GelPro New Life Eco-Pro Anti-Fatigue Mats

1 (20" x 32")	BGPAFM20x32
1 (20" x 48")	BGPAFM20x48
1 (20" x 72")	BGPAFM20x72
1 (24" x 36")	BGPAFM24x36
1 (36" x 60")	BGPAFM36x60



GelPro OR Anti-Fatigue Mats

- Anti-fatigue mats for operating, sterile processing, and surgical spaces, with traction for surgical booties and anti-bacterial and anti-microbial properties

COLOR	QUANTITY	PRODUCT ID#
Columbia Blue	(1) 18" x 24"	GPMM18x24-B
Columbia Blue	(1) 20" x 32"	GPMM20x32-B
Columbia Blue	(1) 20" x 48"	GPMM20x48-B
Columbia Blue	(1) 30" x 72"	GPMM30x72-B
Red	(1) 18" x 24"	GPMM18x24-R
Red	(1) 20" x 32"	GPMM20x32-R
Red	(1) 20" x 48"	GPMM20x48-R
Red	(1) 30" x 72"	GPMM30x72-R



GelPro Medical Disposable Mats

- One-time use, disposable floor mats for added ergonomics in work spaces and step stools

QUANTITY	PRODUCT ID#
1 case (5, 13" x 17" mats)	GPSS-13x17-05
1 case (5, 20" x 30" mats)	GPSS-20X30-05



Fixate™ Anti-Fatigue Mat Gel Pads

- Great to place under mats to prevent slipping

1 pack of 10

AFM-FGP-10



VR™ Sleeve Gloves

- Latex free sleeve gloves provide 29 inches of coverage
- Ultimate protection while cleaning instruments
- Can be reused throughout the day
- Puncture- and tear-resistant



50 SERIES

- Glove length: 28"
- Glove thickness: 11 mil
- Sleeve thickness: 4 mil
- Glove lining: No

SIZE	QUANTITY	PRODUCT ID#
Extra Small	1 box of 50	SGXS-50
Small	1 box of 50	SGS-50
Medium	1 box of 50	SGM-50
Large	1 box of 50	SGL-50
XL	1 box of 50	SGXL-50

80 SERIES

- Glove length: 28"
- Glove thickness: 15 mil
- Sleeve thickness: 4 mil
- Glove lining: Yes

SIZE	QUANTITY	PRODUCT ID#
Extra Small	1 box of 50	SGXS-80
Small	1 box of 50	SGS-80
Medium	1 box of 50	SGM-80
Large	1 box of 50	SGL-80
XL	1 box of 50	SGXL-80



Gloves taking on water?

Strong, durable and latex free protective apparel for splash protection.


PURE
PROCESSING™

NEW

QuickSonic™ Sink System

An ultrasonic and reprocessing sink solution for busy departments

- Height-adjustable to fit each worker's ergonomic needs
- Can be made to fit any sized space and/or instrument processing need
- Performs gross cleaning, flushing, and ultrasonic steps without having the need to transfer the trays from sink to sonic.
- 30"x18" basin dimensions make it perfect for robotics, orthopedic sets, and other larger trays



Ideal for ocular applications:

“

A designated cleaning area and equipment dedicated to the cleaning of intraocular surgical instruments should be identified. Intraocular surgical instruments should be processed separately from general surgical instruments and equipment to reduce the potential for cross-contamination by material or residue from general surgical instruments.

ANSI/AAMI M.2.2.2

”



Healthcare Reprocessing Sinks & Accessories

PureSteel™ SCOPE BUDDY™ PLUS Shelf

QUANTITY	PRODUCT ID#
1	PS-SLF-SBP



Leak tester and Scope Buddy+ shelves to organize cleaning accessories off counters

PureSteel™ Endoscopy Sinks

- Created with flexible & rigid endoscope reprocessing in mind
- Customizable lengths, basin configurations, depths, and more allow sinks to accommodate smaller reprocessing areas
- All the features & accessories offered for sterile processing sinks, but with additional scope reprocessing accessories available:
 - Leak tester shelves
 - Heightened task light magnification
 - Built-in LED basin lighting
 - Solid back walls are easy to clean
 - Temperature gauges to indicate solution temperatures
 - Staging counters and custom staging space solutions
- Standard & custom sizes available



Integrated basin lighting makes leak testing and cleaning easier and more effective



“

At a minimum, two sinks or one sink with two separate basins should be used. Consideration should be given for additional space for delayed processing protocol. Optimally, three sinks or one sink with three separate basins should be used...

ANSI/AAMI ST91: 2021 4.3.2 Sinks and accessories

”




PureSteel™ Healthcare Reprocessing Sinks

- Customizable, complete-solution sink engineered around instrument reprocessing
- Electric push button allows precise, personal height adjustment; two buttons on each sink
- Ergonomic, slanted wrist edges on edge of sink basin remove pinch points during reprocessing
- Optional posterior back wall to organize cleaning accessories; adjusts up & down as the sink adjusts height
- Open configuration storage shelf underneath the sink raises and lowers to allow mopping and cleaning, and raise items for staff
- Inline sonic systems available for added ultrasonic capacity
- Adheres to industry guidelines for reprocessing and ergonomics
- Standard sizes: 96", 108", or 120"
- Custom sizes available

Slanted Ergonomic Wrist Rests

QUANTITY	PRODUCT ID#
Per Basin	SNK-ERG-SLANT



Remove painful pinch points during cleaning

PureSteel Xi Reprocessing Sink

- Fixed or height-adjustable
- Customizable packages available
- Allows departments to create a dedicated reprocessing area
- Small footprint
- Can accommodate robotics, MIS, and larger scopes



Basin sizes optimized for larger instrumentation, such as robotics and orthopedic



QUANTITY	PRODUCT ID#
QTY 1 (height-adjustable)	PSXI-HA-SA
QTY 1 (fixed height)	PSXI-SA



PureSteel™ EZ Gate Hinge™ Pegboards



- Optional, hinge pegboards help technicians and nurses clean behind pegboards or back walls
- In tight alcoves or installations, hinge pegboards provide control to move accessories, access hardware, utilities, and more



PureSteel™ Inline Sonic Sink System

Great for adding capacity for non-lumened devices in small-department footprints.

- Integrated ultrasonic system in any basin of choice for non-complex sonic cycles
- Customizable sonic cycles & times
- Removable for easy servicing

QUANTITY	PRODUCT ID#
30" Basin	PS-INLINE-SONIC-301810B



PureSteel™ LED Basin Lights

- In-basin lighting aids in spotting bubbles during leak testing
- Enhances lighting at the sink, making it easier to perform cleaning tasks and check for quality

QUANTITY	PRODUCT ID#
2 lights	PS-LED-BSN-02
3 lights	PS-LED-BSN-03



Sink Staging Counter Storage Shelf

- Built-in storage shelves keep staging counter organized and easily accessible

1	PS-SC-STSH-C
---	--------------



PureSteel™ Temperature Gauge

- Integrated temperature gauge that provides consistent, constant temperature readings for PureSteel™ Sinks

1	PS-THERMO-GAUGE
---	-----------------




Paladin® Evolution Rail System for Sinks

- Mounts directly to sink pegboards allowing users to easily 'snap', move and re-organize accessories to different locations, or replace accessories
- Great for EZ Hinge Pegboards, or sinks hard mounted to the floor or within tight areas



QT 2, 2ft rail kit	PS-SK-EVOLUTION-2FTKIT
--------------------	------------------------



PureSteel™ Splash Guard

- Provides a clear shield to cover a sink basin, reducing a user's contact with aerosols and contaminated reprocessing fluids. Splash guards can be lifted up when not in use.

1	PS-SGP-01
---	-----------



PureSteel™ Staging Panel

- Create optional and flexible arrangements of additional counter space when staging room is required.
- Each Staging Panel can be placed onto a sink basin to create extended work surfaces for reprocessing.

1	PS-SC-01
---	----------

NEW



PureSteel™ Copper Clean Anti-Microbial Wrist Rest

- Provides continuous anti-microbial, self-sanitization properties in high-touch environments.
- Copper can resist corrosion and eliminate up to 99.9% of six common microorganisms in less than two hours, including: MRSA, VRE, C. Difficile, E.coli, Influenza A and Norovirus

QUANTITY	PRODUCT ID#
1	PS-COPPER-SLANT-28B
1	PS-COPPER-SLANT-30B


NEW



PureSteel Workflow Plate

- Plates indicate location for reprocessing steps: soak, wash, rinse, final rinse, leak, flush or test options.
- Removable, cleanable alternative to sink basin stickers

1, Wash	PS-C-WS-ETCH
1, Soak	PS-C-SK-ETCH
1, Rinse	PS-C-RS-ETCH
1, Leak	PS-C-LK-ETCH
1, Flush	PS-C-FL-ETCH
1, Test	PS-C-TS-ETCH
1, Ocular Sink	PS-C-EYE-ETCH




PureSteel Volume Indicator Plate

- Removable, cleanable alternative to sink basin stickers
- Custom made to any stainless steel sink basin size



1	PS-C-VL-ETCH
---	--------------



PureSteel™ Stainless Steel Solution Tank

- Provide a separate, dedicated space for chemicals and solutions when reprocessing

1, right-side tank	PS-SSST-16BM-RIGHT
1, left-side tank	PS-SSST-16BM-LEFT
1, Xi Sink size	PS-SSST-16m-XISINK



PureSteel Rinse Sprayer

- DI rinse sprayer provides consistent and copious rinsing of instrument sets and trays using treated water.

QUANTITY	PRODUCT ID#
1 rinse sprayer with tubing and mounting hook	ASG-11
1 rinse sprayer with adapter base and mounting hook	ASG-11-LT-V250-T1
1 all-in-one spigot and sprayer	PS-DI-SPIGOTGUN



PureSteel Telescoping Drains

- Telescoping drains are made of chrome plated brass construction with smooth walls
- Creates a durable, long-lasting drain system for vertically adjusting equipment

QUANTITY	PRODUCT ID#
1 drain	T1512



Pure LandingPad Silicone Liner

- 20" x 20" silicone sink landing pad to protect delicate instruments, such as endoscopes
- For use in sink basins or counter tops



1 liner	PLP20x20-01
---------	-------------



Pure LandingPad™ Surface Liner

- Protects delicate devices and cushions surfaces where expensive endoscopes and other devices are staged for cleaning
- A 3/4-inch thick nitrile rubber foam blend resists hospital-grade cleaning chemicals and disinfectants and will not absorb liquid
- 18" x 24"



1 pack (3 liners)	PLPSL-1824-03
-------------------	---------------



PureSteel Light Hood

- Durable light hood provides high-quality lighting at the sink for optimal visualization of surfaces and soils.

1 (4' sink hood)	T-SK-PB-HOOD-04
1 (2' sink hood)	T-SK-PB-HOOD-02
1 (4' table hood)	T-TB-PB-HOOD-04
1 (2' table hood)	T-TB-PB-HOOD-02



PureSteel™ UV Blacklight LED Light Hood

- UV blacklights to aid in hands-free detection of residual cement from orthopedic devices
- Can be used independently of regular overhead LED lighting




1 (2' sink hood)	PS-UVLEDBLACK-02
------------------	------------------



Brush Bath and Shelf

- Allows instruments to soak on a sink pegboard or back wall with a built-in draining system
- Made from stainless steel

QUANTITY	PRODUCT ID#
1 (12")	BBATH-12



PureSteel Brush Holder and Organizer

- Organizes cleaning brushes at an easy to reach, ergonomic height for users
- Integrates with existing pegboards and allows up to 8 brushes to be hung

1	PS-BRRK
---	---------



PureSoak™ Timer

- Water-resistant timer for use on sinks and soak stations
- Highly visible LED digits
- Programmable minutes and seconds which count-down; alarm activates when the programmed time has elapsed



1	PSWPT
---	-------



Automated Endoscope Reprocessing Rack

- Hangs on a wall; can plug adapter in close to the automated endoscope reprocessor to keep tubes organized
- Can accommodate up to 11 AER connectors

1 (Steris, Olympus, Evotech)	AER-TSO-01
1 (Medivators)	AER-TSO-02



Bluetooth Speaker

- Connects by bluetooth to phones, computers
- Small profile
- Great for reprocessing sinks, tables

1	PS-BTSPEAKER
---	--------------



PureSteel Pegboard Panels


- Individual pegboard panels for back wall organization
- Install accessories at a comfortable working height and keep counters organized



1 (2' x 2' panel)	SSS-PB-2x2
1 (2' x 4' panel)	SSS-PB-2x4

Leak Tester Shelves

- Shelves which allow for storage of leak testers onto back walls and pegboards



QUANTITY	PRODUCT ID#
1, Small Size	PS-SM-LTS
1, Olympus	P S-SLF-OKV6



PureSteel Glove Box Holder

1 (12" L x 6" H x 4" D)	GBH-01
1 (16" L x 5" H x 2" D)	GBH-02



PureSteel Wire Baskets

- Organize reprocessing equipment and supplies onto a pegboard at ergonomic height for staff, and create more counter space

1 (12"X12"X8")	HP6203
1 (24"X10"X5")	HD26298
1 (12"X6"X6")	HP6209
1 (15"X12"X5/3")	HD26302
1 (24"X12"X4")	HD26297
1 (12"X12"X4")	HP6205
1 (12"X12"X8"/4")	HD26300



PureSteel 32oz Bottle Holder

- Conveniently keeps bottles organized on a sink pegboard or back wall
- Holds one or two 32oz bottles
- Made from stainless steel

1 (single holder)	SBH-01
1 (double holder)	DBH-02



PureSteel™ Gallon Holder

- Conveniently keeps gallon bottles organized on a sink pegboard or back wall
- Holds one or two 1-gallon bottles
- Made from stainless steel

1 (single holder)	SGH-01
1 (double holder)	DGH-01



Mobile Solutions & Carts

PureSteel Mobile Soaking Station

- Create dedicated spaces for soaking during peak reprocessing time
- Built on casters for easy transport
- Create extra sink spaces for reprocessing, or establish dedicated areas for specific instruments
- Easily accommodates robotic instruments, scopes, or other longer instrumentation
- Back wall pegboard organizes equipment at the area of manual cleaning
- Features ball valve for draining to floor drain; sink hoppers



Meet robotic soaking IFU




QUANTITY	PRODUCT ID#
1 (34"x22")	PS-MSSPB-34x22
1 (optional faucet kit)	FSNK-01
1 (replacement tubing)	MSS-TUBING
1 (power strip)	MSS-PS-01



PureStation Mobile Cart

- Enables point-of-use cleaning
- Can be stored out of the way and brought out at peak processing hours as an additional sink station
- Can be used to establish a separate dedicated area for specified instruments, such as eyes
- *Must be used with the PureStation Sink Insert*



Meet ocular guidelines



1	PSCART-1626
1 (ball valve kit for draining)	PSBV-1KIT



PureStation Mobile Cart Revital-Ox® RESERT® Bracket

- Adapter for Mobile Cart frame to fit Revital-Ox RESERT HLD soaking containers.
- Provides a mobile platform for HLD soaking.
- Frees up counter space by providing an additional area for manual bins
- **Must be used with PureStation Mobile Cart (PSCART-1626)**



1 bracket	PSCART-1626-RESERT
-----------	--------------------

PureSteel XL Transport Container Cart

- Provides a secure, dedicated space to horizontally soak instruments
- Also serves as a sealed soaking station for facilities lacking adequate soaking space
- Organizational back wall allows for installation of the FlexiPump Independent Flushing System and other accessories
- Under-shelf allows for additional storage
- XL Transport Containers purchased separately



Meet
robotic
soaking
IFU



QUANTITY	PRODUCT ID#
1 Cart	XL-TRSP-CART



PureSteel XL Transport Container Tiered Cart

- Provides a secure, dedicated space to horizontally soak or transport in up to four InstruSafe XL Transport Containers
- XL Transport Containers purchased separately



1 Cart	XL-4TIER-CART
--------	---------------



InstruSafe® XL Transport Container

- Creates a dedicated space for horizontally soaking robotic instruments
- Dedicated dry erase "soak start time" space on cover
- Inner gasket and strong silicone latches reduce spills
- Quick-drain feature — no need to tip container to drain



1 (w/ drain)	XL-TRSP-DR
1 (w/out drain)	XL-TRSP



Prep and Pack Tables & Accessories


PureSteel Ergonomic WorkStations

- Optional tabletop light provides illumination for inspection of sterile wrap materials
- Height-adjustability allows technicians to perform tasks while standing or seated
- Wide catalog of accessories, from monitors to shelving, to tailor tables to their needs and use
- Optional, fully customizable back walls organize accessories, power cords and computers at the most ergonomic height for technicians; adjust up & down as the table height adjusts
- Wide counter space provides more room for larger instruments and wraps
- Equipped with eight electrical outlets and cord management
- Standard sizes: 60", 72" and 78" lengths
- Custom sizes available




PureSteel™ Hughes Seal & Stow Heat Sealer Shelf

QTY	PRODUCT ID#
1	PS-TB-HEATSEALSLF



Stores large, heavy seat sealers off of countertops.



PureSteel EZ Slide Scale

- Pull-out scale that allows technicians to weigh trays, instruments, and items during packaging and assembly
- Helps prevent over-packing or exceeding weight limits, improves technician ergonomics, helps ensure compliance.
- Scale dimensions: 10"L x 10"D
- Display dimensions: 9"L x 8.5"H x 3"D

QUANTITY	PRODUCT ID#
1	PS-ERG TB-EZ-SCALE
1 scale, standalone	PS-SCALE




Table Acrylic Kit w/ Ruler

QTY 1 Acrylic Kit

PRODUCT ID#

PS-TB-LGHT-ACRYLICINSERT-RL




Ergonomic Table Pull Out Shelf

QTY

PRODUCT ID#

1

PS-TB-PSHLF-12W



PureSteel Signature Series Ergonomic WorkStations

- For departments looking to add the same ergonomics and productivity benefits with fewer customization options
- Available in 60"x32" and 70"x32" dimensions, height-adjustable, with 24"H pegboards
- Same library of organization-enhancing accessories and tools



PureSteel Wrap Inspection Table

- Rack allows for organization of multiple-sized wraps
- Tabletop light assists in visual inspection of wrap materials
- Optional height-adjustability relieves strain on technicians and meets ergonomic guidelines
- Custom sizes available

Customizable options and back walls available



PureSteel Mini Inspection Table

- Tabletop light assists in visual inspection of wrap materials
- Small footprint
- Casters allow for mobility
- Bottom shelf creates space for storage of materials



NEW

PureSteel Mobile Scale Cart

- Small, mobile scale which can be easily moved throughout a department for weighing trays, instruments, and more.
- Operates on a rechargeable battery, allowing the cart to operate anywhere in the department, regardless of access to an outlet.
- Cart outer dimensions: 24"L x 24"D
- Scale dimensions: 10"L x 10"D
- Display dimensions: 9"L x 8.5"H x 3"D



1 cart PS-MS-CART



PureSteel Prep and Pack Organization Rail

- Hang peel pouch rolls, stickers, and other sterilization materials on a rail for easy dispensing and filling

QUANTITY	PRODUCT ID#
1 rail shelf	PS-TB-RAIL



PureSteel Clear Shelf Bins with Lid


- Clear bins allow users to easily identify bin contents
- Featured lids prevent dust & lint from accumulating on contents
- Compatible with our existing rail system

QUANTITY	PRODUCT ID#
1 large bin	PS-PSB-4X10.87X5
1 standard bin	PS-PSB-4X7.5X3

PureSteel Shelf Bins and Plastic Shelf Bin Rail

- Organization bins for prep and pack tables or pegboard panels
- Requires rail for mounting
- Available in: 5.5" W x 11" L x 5" H

1, 24" blue bin rail	PS-PSBR-24X3
1, 5.5"W bin	PS-PSB-5.5X11X5



PureSteel Peel Pouch Rack Shelf

- Vertically stores and organizes pre-cut peel pouches on a back wall shelf

1, 5 slot peel pouch	PS-PPRS-05
1, 8 slot peel pouch	PS-PPRS-08



PureSteel Lighted Prep and Pack Shelves

- Available in 2', 3' or 4'
- All shelves include built-in light bars for increased illumination underneath


1 (24" shelf)	PS-SLF-LGHT-24X12
1 (36" shelf)	PS-SLF-LGHT-36X12
1 (48" shelf)	PS-SLF-LGHT-48X12



PureSteel Pegboard Hooks

- Hang accessories from pegboards; hardware-free installation for rearrangeable layouts


Pack of 5	PS-HOOK
-----------	---------



Tabletop Ruler

- Accurately measure packaging and assembly materials
- Available in custom lengths
- Can be printed on acrylic replacement panels

1	PS-TB-RULER
---	-------------



PureSteel Tape Dispenser

- Stores and protects sterilization tape from exposure to lighting

1	PS-TAPEDS
---	-----------



PureSteel Peel Pouch Box Holder

- Allows departments to create vertical storage of peel pouch boxes on the sides of tables for convenient pulling, refilling and storage.

32" L Peel Pouch Box Holder	PS-PPH-32
Rail Kit	PS-TB-EVOLUTION-2FT

PureSteel Specialty Impact Tables


The operating theatre is as busy as ever! Surgical cases are at an all-time high, specialty devices are diversifying workflows and complicating compliance requirements while sterile processing is ramping up productivity with strained resources.

- PureSteel Specialty Impact Tables are designed to help departments by:
- Delivering productivity by focusing technicians and tasks in centralized areas and focused equipment set-ups
- Ensuring compliance with unique workflows, IFU requirements and new and updating standards
- Organizing limited spaces, and keeping areas free of clutter, and protecting devices and tools



PureSteel™ Robotics Workstation

- Boost productivity and compliance with specialty devices
- The PureSteel Robotics WorkStation dedicates space and resources for robotic packaging areas.
- Features:
 - Hughes Seal & Stow Heat Sealer Shelf frees up countertop space
 - Side storage for robotic-sized peel pouches at technicians' fingertips
 - Pegboard backwall integrates 4.75X task lighting, light hoods, lighted shelving for printers, tracking systems, and other key tools.



PureSteel™ Low-Temp Sterilization Workstation

- Organize and ensure compliance across different sterilization methods
- For departments seeking to add organization, efficiency and compliance near low-temperature sterilization processes.
- Features:
 - Integrated scale for weighing trays
 - Back pegboard hosts a variety of accessories, including:
 - Ergotron mounts
 - Shelves for printers
 - Basket bins for indicators, locks, challenge packs, etc.
 - Tape dispenser to protect sterilization tape from over exposure
 - Under-table storage for mounting computers, printers, and more!




PureSteel™ Loaner WorkStation

- Organize and set standards with loaner delivery and pick-up processes
- For departments with high volumes of loaner trays, or when loaner tray processing requires more oversight to help increase control over a department's loaner program.
 - Integrated scale to ensure trays are within guidelines for tray weight
 - Pegboard to mount and organize critical accessories, such as:
 - Cameras
 - Tracking system computers
 - Scanners
 - Task lights dedicated for detailed instrument inspection



**Quality Assurance
Solutions**



PureClear™ Inspection Scopes

- Allows technicians and nurses to visually inspect lumens and devices for debris and damage
- Inspects channels, connection points, distal tips, and spaces behind elevator channels
- USB Control Box provides video and photo capture for documentation
- HDMI Control Box gives inspection capabilities without the need of a PC or computer

QUANTITY	PRODUCT ID#
BORESCOPIES	
(1) 1.83mm, Length 60cm	CB2-60
(1) 1.83mm, Length 110cm	CB2-110
(1) 1.83mm, Length 200cm	CB2-200
(1) 1.06mm, Length 110cm	CB1-110
CONTROL MODULES	
PureClear™ Dual Inspection Control Module, USB Output	CD-USB
PureClear™ Dual Inspection Control Module, HDMI Output	CD-HDMI
MOUNTS	
1 (mount)	CD-FM

PureClear Visualization Station

- Provide a dedicated space for visual inspection or work
- Can be height-adjustable for greater staff ergonomics
- Can be used on the decontamination side of the department to visually inspect items before continuing with further reprocessing
- Easily integrate monitors for larger visual displays
- Standard visualization station available in 36" x 30"
- Mini visualization station is available in 24" x 24"



QUANTITY	PRODUCT ID#
1 (fixed height)	PC-VSSTATION
1 (height-adjustable)	PC-HA-VSSTATION
1 (mini fixed height)	PC-MINI-VSSTATION
1 (mini height-adjustable)	PC-HA-MINI-VSSTATION



Intel NUC Mini PC Kit

- Kit provides PC hard drive capabilities for touch pads and monitors.

QUANTITY	PRODUCT ID#
1	PC-NUC-PC-KIT



PureSteel Magnification Task Light

- Up to 3x magnification
- Features on/off switch and dimmable lighting
- 31" or 39" arm lengths
- Lights come with c-clamps standard; L-bracket option sold separately
- Optional, 6 diopter lens increases total task light visibility to 2.25x magnification



1 (39"L)	PS-LED-MG-LGHT-SS-39
1 (31"L)	PS-LED-MG-LGHT-SS-31
1, L-Bracket	PS-LED-MG-WB
1 (6D lens)	PS-LED-MG-LENS-6D



PureSteel 4.75x Magnification Task Light

- Meets 4x magnification requirements
- 2 removable, swappable lenses (no tools required)
- LED bulb; est. 20,000 hour lifespan
- Touch-sensitive brightness controls
- 2.25x lenses also included



1	PS-LGHT-475X
1, replacement 4.75x len	PS-LGHT-475X-LENS



PureSteel Halogen Magnifier Task Light

- 2x magnifier assists staff in visual inspection of instrument tips and surfaces
- Targeted illumination with a 27" flexible gooseneck arm
- Water resistant to work in moist environments

FOR PEGBOARD INSTALLATION	
1	PS-HAL-LIGHT
MOUNTING OPTIONS (non-pegboard)	
1, L-Bracket	PS-LGHT-BRCKT
1, C-Clamp	C-CLAMP-LGHT
REPLACEMENT PIECES	
1, magnifier lens	PS-HAL-MGF



All-in-One Ergotron® Mount Arm

- Mounts a monitor, scanner holder, and keyboard and mouse platform for IT stations in one arm

QTY	PRODUCT ID#
1	PC-VS-ARM-S-KYBD



Ergotron® Wall Mount Arm

- Mounts a screen directly to a wall; only mounts a monitor.

1	PC-VS-ARM
---	-----------



Ergotron® LX Wall Mount Keyboard Arm

- Mounts keyboards and a mouse pad to a wall, creating space on a work surface. Does not mount a screen.

1	PC-VS-KYBD
---	------------



Ergotron Lockable Tablet Mount

- Converts a monitor arm to a tablet mount for smaller screens or touch tablets.

1	PC-VS-TABLET
---	--------------



Peerless Tilt Wall Mount

- Mounts a screen directly to a wall. An arm-less alternative to monitor or screen mounts.

1	PC-VS-ST630P
---	--------------



Ergotron® Wall Mounted Handheld Scanner Mount

- Mounts scanner holders directly to a wall

1	PC-WALL-HSM
---	-------------



Ergotron® VESA Scanner Shelf

- Mounts scanner shelves to the side of monitor mounts, taking scanners off of work surfaces

QUANTITY

PRODUCT ID#

1

PC-VS-VSS



Ergotron® Dual Monitor & Handle Kit

- Adapts a single monitor arm mount to attach a second screen. Must be purchased with another monitor arm mount.

1

PC-VS-DUAL



Ergotron® Mini Desk Arm

- Mounts monitor arms to the sides of a surface or desk for when back walls are unavailable

1

PC-VS-ARM-DESK



Monitor Pole Arm

- Allows tablets and screens to be mounted on a pole

1

PC-VS-ARM-POLE



Elo® Touchscreen Monitor

- Well-suited for use in high traffic or harsh environments

1 (15" monitor)

PS-VS-ELO-15

1 (20" monitor)

PS-VS-ELO-20

1 (24" monitor)

PS-VS-ELO-24



PURE
PROCESSING™

2200 Ogden Ave, Suite 100, Lisle, Illinois 60532
877-718-6868 | pure-processing.com

Raise the Standard